



We Help The World Play Safe



BRAND NAME: CUPID

PRODUCT FEATURES

Full product name and branding		
Description	Plain, Dotted, Ribbed, Multi-textured, Super-dotted, Super-dotted, Extra Time, Xtra Large, Xtra Strong	
Length	≥170mm, ≥180mm, ≥190mm	
Width	49+/-2mm, 53+/-2mm, 56+/-2mm	
Shape (parallel, anatomical, contoured, other)	Parallel, Teat-ended	
If other, please specify	As per requirement	
Colour	Available In Different Colors (*As Per Customer Requirement)	
Flavour	Available In Different Flavours (*As Per Customer Requirement)	
Lubricant (type and quantity per condom) Silicone Oil (Can be adjust		



as per Customer requirement)

SQUARE PACKAGING		
Type of packaging	3Pcs, 6pcs , 10pcs, 12pcs,	
Instructions for use / leaflet	1 IFU or Backside Printing	
Dimensions	As per customer requirement	
Paper Board Specs	As per customer requirement	
Celophane / shrinkwrap	As per customer requirement	
Artwork No. Of Colors	CMYK, SPECIAL	
Artwork Format	Ai/CDR	

RECTANGULAR PACKAGING		
Type of packaging	3Pcs,	
Instructions for use / leaflet	Backside Printing	
Dimensions	As per customer requirement	
Paper Board Specs	As per customer requirement	
Celophane / shrinkwrap	As per customer requirement	
Artwork No. Of Colors	CMYK, SPECIAL	
Artwork Format	Ai/CDR	









CIRCULAR PACKAGING		BULK PACKAGING	
Type of packaging	1Pcs, 3Pcs, 6Pcs,10Pcs, 12Pcs	Type of packaging	144Pcs
Instructions for use / leaflet	1 IFU or Backside Printing	Instructions for use / leaflet	1 to 10 IFU
Dimensions	As per customer requirement	Dimensions	As per customer requirement
Paper Board Specs	As per customer requirement	Paper Board Specs	As per customer requirement
Celophane / shrinkwrap	As per customer requirement	Celophane / shrinkwrap	As per customer requirement
Artwork No. Of Colors	CMYK, SPECIAL	Artwork No. Of Colors	CMYK, SPECIAL
Artwork Format	Ai/CDR	Artwork Format	Ai/CDR





IT'S TIME TO FEEL THE CHANGE



BRAND NAME: CUPID ANGEL, D'SIRE

PRODUCT FEATURES

Description	Version-1, Version-2,
Length	145.0-165.0, 125.0-145.0mm
Nominal Width	75+/-2mm
Shape	OCTAGON
If other, please specify	As per requirement
Colour	Available In Different Colors (*As Per Customer Requirement)
Coloui	Available In Different Flavours (*As Per Customer Requirement)





- Unique by its design.
- Female Condom available in Various Colors.
- Female Condom available with Vanilla Scent.
- Economical / Cost Effective & Easy to use.
- Made from finest Quality Latex.

Female condom is women controlled dual protection device to prevent the spread of STI's including HIV and to prevent unintended pregnancies.

TYPES OF PACKAGING

PACKAGED		
Type of packaging	1Pcs, 2pcs	
Instructions for use / leaflet	1 IFU	
Dimensions	As per customer requirement	
Paper Board Specs	As per customer requirement	
Celophane / shrinkwrap	As per customer requirement	
Artwork No. Of Colors	CMYK, SPECIAL	
Artwork Format	Ai/CDR	

BULK		
Type of packaging	100Pcs	
Instructions for use / leaflet	10 IFU	
Dimensions	As per customer requirement	
Paper Board Specs	As per customer requirement	
Celophane / shrinkwrap	As per customer requirement	
Artwork No. Of Colors	CMYK, SPECIAL	
Artwork Format	Ai/CDR	



BRAND NAME: CUPID LUBRICANT

Full product name and branding	WATER BASED LUBRICANT	
Description	WATER BASED LUBRICANT	
Osmolality	Less than or equal to 1200 mOSm/ kg	
Glycerol Content	Less than or equal to 9.9 % mass fraction	
Propylene Content	Less than or equal to 8.3 mass fraction, pH 5.5-7.0, Viscosity: 8000 - 20000 cps	
INGREDIENTS:	Carbopol Glycerin Sodium Benzoate Carbopol NaOH Purified Water	
If other, please specify	As per requirement	
PROPERTIES	Water soluble Viscous Non Greasy Clear Transparent Flavored / Non Flavored	
Flavour	As per customer requirement	

SACHETS: 4ML, 5ML, 10ML,		
Type of packaging	Packaged & Bulk	
Instructions for use / leaflet	1 IFU	
Dimensions	As per customer requirement	
Paper Board Specs	As per customer requirement	
Celophane / shrinkwrap	NA	
Artwork No. Of Colors	CMYK, SPECIAL	
Artwork Format	Ai/CDR	

WATER **BASED** LUBRICANT

TUBES: 20ML, 30ML, 50ML, 82ML, 114ML	
Type of packaging	Packaged & Bulk
Instructions for use / leaflet	1 IFU
Dimensions	As per customer requirement
Paper Board Specs	As per customer requirement
Celophane / shrinkwrap	NA
Artwork No. Of Colors	CMYK, SPECIAL
Artwork Format	Ai/CDR

Lubricant Applications

Provides very effective lubrication when women's are suffering vaginal dryness. Stimulate the genitals during sexual intercourse. It does not contain any kind of spermicidal content. Safer to use & does not have any side effects.

ADVANTAGES OF LUBRICANT:

- Unique in its formulation.
- Available in various flavors.
- Spreads readily, adheres well, Harmless to human tissue.
- Does not damage rubber appliances.
- User friendly pack sizes.
- Economical/ Cost Effective. Made from finest quality Raw Materials meeting standard specifications. Compatible with natural latex condoms, polyurethane and poly-isoprene condoms.

SPECIFICATION

- pH (at 25 °C): 5.5 to 7.0
- Osmolality (mOsm/kg): Not More Than 1200
- Viscosity (Cps.) At 25 °C: 8000 to 20000
- Bio-Burden Level: Below 100 cfu per gram & Pathogens are absent

• Shelf Life: 3 years from the date of manufacturing













We Help The World Play Safe

Factory & Registered Office:

Cupid Limited, A-68, M.I.D.C.(Malegaon) Sinnar, Nashik-422113, Maharashtra, INDIA. Tel:+91-2551-230280 / 230772

info@cupidlimited.com | www.cupidlimited.com | www.cupidangel.in