



**CUPID**  
LIMITED

*We Help The World Play Safe*



**LUBRICANT**  
**WATER BASED**

Brochure



<b>Full product name and branding</b>	<b>WATER BASED LUBRICANT</b>
Description	<b>WATER BASED LUBRICANT</b>
Osmolality	Less than or equal to 1200 mOSm/ kg
Glycerol Content	Less than or equal to 9.9 % mass fraction
Propylene Content	Less than or equal to 8.3 mass fraction, pH 5.5-7.0, Viscosity: 8000 - 20000 cps
<b>INGREDIENTS:</b>	<ul style="list-style-type: none"> <li>• Carbopol</li> <li>• Glycerin</li> <li>• Sodium Benzoate</li> <li>• EDTA</li> <li>• NaOH</li> <li>• Purified Water</li> </ul>
If other, please specify	<b>As per requirement</b>
<b>PROPERTIES</b>	<ul style="list-style-type: none"> <li>• Water soluble</li> <li>• Viscous</li> <li>• Clear Transparent</li> <li>• Non Irritant</li> <li>• Non Greasy</li> <li>• Non staining</li> <li>• Flavored / Non Flavored</li> </ul>
Flavour	<b>As per customer requirement</b>

### SACHETS : 4ML, 5ML, 10ML,

Type of packaging	<b>Packaged &amp; Bulk</b>
Instructions for use / leaflet	<b>1 IFU</b>
Dimensions	As per customer requirement
Paper Board Specs	As per customer requirement
Cellophane / shrinkwrap	NA
Artwork No. Of Colors	CMYK, SPECIAL
Artwork Format	Ai/CDR

### TUBES : 20ML, 30ML, 50ML, 82ML, 114ML

Type of packaging	<b>Packaged &amp; Bulk</b>
Instructions for use / leaflet	<b>1 IFU</b>
Dimensions	As per customer requirement
Paper Board Specs	As per customer requirement
Cellophane / shrinkwrap	NA
Artwork No. Of Colors	CMYK, SPECIAL
Artwork Format	Ai/CDR

### ADVANTAGES OF LUBRICANT :

- Unique in its formulation.
  - Available in various flavors.
  - Spreads readily, adheres well, Harmless to human tissue.
  - Does not damage rubber appliances.
  - User friendly pack sizes.
  - Economical/ Cost Effective.
- Made from finest quality Raw Materials meeting standard specifications. Compatible with natural latex condoms, polyurethane and poly-isoprene condoms.

### SPECIFICATION

- pH (at 25 °C): 5.5 to 7.0
- Osmolality (mOsm/kg) : Not More Than 1200
- Viscosity (Cps.) At 25 °C : 8000 to 20000
- Bio-Burden Level: Below 100 cfu per gram & Pathogens are absent
- Shelf Life: 3 years from the date of manufacturing

### Lubricant Applications

Provides very effective lubrication when women's are suffering vaginal dryness. Stimulate the genitals during sexual intercourse. It does not contain any kind of spermicidal content. Safer to use & does not have any side effects.



Follow us:    



**CUPID**  
LIMITED

---

## **Cupid Limited**

Factory Address : A-68, M.I.D.C. (Malegaon) Sinnar,  
Nashik-422 113, Maharashtra, INDIA.

Tel : +91-2551-230280 / 230772

E : [info@cupidlimited.com](mailto:info@cupidlimited.com)

W : [www.cupidlimited.com](http://www.cupidlimited.com)